

DeepChek® Assays: Your best ally for HIV genotyping & Drug Resistance Determination

From Capillary Electrophoresis to Deep Sequencing: An Improved HIV-1 Drug Resistance Assessment Solution Using In Vitro Diagnostic (IVD) Assays and Software

Sofiane Mohamed (1,*), Ronan Boulmé (1) and Chalom Sayada (2)

1 ABL Diagnostics, 57140 Woippy, France
2 Advanced Biological Laboratories (ABL), 2550 Luxembourg, Luxembourg
* Correspondence: s.mohamed@ablsa.comcass

Objective:

The objective of this study was to evaluate the performances of CE-IVD assays for HIV-1 drug-resistance assessment both for **target-specific** and **whole-genome sequencing**, using standardized end-to-end solution platforms.

Method:




Figure 1

Figure 2

Localization of the CE and NGS primers for HIV drug resistance:

- In blue, CE primers for reverse-transcriptase, protease, and integrase regions.
- In purple, NGS primers for the whole-genome HIV (Snappene Software Version 5.25.5)




Table 1

A comparison of the NGS protocol and CE sequencing is summarized in Table 1. To process 24 samples (RT, PR, and INT) by CE sequencing, the time required for waiting, sample preparation, and overall time to result were 81.0 h, 2.0 h, and 83.0 h, respectively. The full workflow of database processing, analysis, and reporting using DeepChek®-HIV was more than 20 min per sample. To process 24 samples (RT, PR, and INT) by NGS, the time required for waiting, sample preparation, and overall time to result were 27.0 h, 4 h, and 31.0 h, respectively. The full workflow of database processing, analysis and reporting using DeepChek®-HIV was less than 2 min per sample.

Comparison of NGS and CE methods for the HIV-1 target-specific sequencing.

Steps	NGS	Time/24 Samples (h)	CE	Time/24 Samples (h)
Sample preparation	RNA extraction kit	1.0	RNA extraction kit	1.0
Amplification	RT-PCR	4	RT-PCR	4
Purification	Beads Purification	0.75	Enzymatic purification	0.2
Quantitation	Quality control (TapeStation)	0.2	-	-
	Normalization (Qubit)	0.5	-	-
Library/sequencing reaction	Library preparation	4	Sequencing reaction	2.5
Dilution	Dilution and pooling	20	Sequencing with SeqStudio	72
Sequencing	Sequencing		4-capillary	
	FastQ files		ABI files	
Data analysis	DeepChek® using ANRS, HIVdb, etc.	0.2	DeepChek® using ANRS, HIVdb, etc.	1.0
	Handling time	4	Handling time	2
Result	Waiting time	27	Waiting time	81
	Time to result	31	Time to result	83
Price	Reagent cost \$/sample	100-150 *	Reagent cost/sample	80
Sensitivity		1 to 3%		20%

* Including extraction, PCR, library preparation, indexes, sequencing, and software.

Conclusions:

The use of whole-genome sequencing is an additional and complementary tool to detect mutations in newly infected untreated patients and heavily experienced patients, both with higher HIV-1 viral-load profiles, to offer new insight and treatment strategies, especially using the new HIV-1 capsid/maturation inhibitors and to assess the potential clinical impact of mutations in the HIV-1 genome outside of the usual HIV-1 targets (RT/PR and INT).



ABL Diagnostics S.A.

72c, rue de Thionville – 57140 WOIPPY, France

contact@abldiagnostics.com

<https://www.abldiagnostics.com>

Phone :+334 9105 1878